



DOBSON ORTHO LABORATORY, INC.
 2130 Hilton Drive, Suite B
 Gainesville, GA 30501
 Tel: (888) 599-7890 GA: (470) 290-8601
 Fax: (470) 290-8600
dobsonortho.com

Sintered / Printed Metal Rx

LAB USE ONLY

EXPANSION APPLIANCES Upper Lower

Hyrax RPE Standard (4-arm)
Dentaurum Palex Neo

Hyrax RPE Ratchet (4-arm)
Forestadent Snap Lock

Mini RPE Standard (2-arm)
Dentaurum Variety SP

Add Habit Crib to Design - Select Type
 With Spikes No Spikes

HOLDING APPLIANCES Upper Lower

Transpalatal Arch - no loop

Lingual Arch - no loops

Space Maintainer _____

HABIT APPLIANCES - Include counter model

Open Spike Design
 Closed Cage Design

ACCESSORIES

Face Mask Hooks - Select Type
 Vertical / Intrusion Horizontal / Protraction

Rests on 4/d's 5/e's

Extensions to 6's/ 7's

METAL ALLOY NOTICE
 All sintered parts are printed using a chromium cobalt alloy. Expansion screws and some accessories may still be stainless steel. Please review the *Printed Metal Delivery Protocol* on *DobsonOrtho.com* BEFORE you seat your first appliance as the bond to chrome may be stronger than stainless steel.

LEAD TIME
 Each appliance is 3D designed using specialized CAD software and then are 3D printed using a sintering printer. This process can take up to 8 days before final welding and/or acrylic. We ask that you please account for a little extra time for this product than you would expect for our standard fixed metal appliances.

ADJUSTMENT WARNING NOTICE
 Due to the rigidity of the chromium cobalt alloy, we do not put any adjustment loops in our Sintered Metal Appliances as trying to adjust any of the sintered wires WILL break the appliance and void any warranty.

LIMITED WARRANTY
 Due to the high cost of the sintering process, we reserve the right to evaluate any models or scans when a remake has been requested. PLEASE VERIFY YOUR SCANS ARE COMPLETE BEFORE SUBMITTING FILES TO OUR LAB.

DOCTOR _____ ACCT# _____

ADDRESS _____

CITY _____ STATE _____ ZIP _____

PHONE _____ FAX _____

EMAIL _____

PATIENT NAME _____

DATE SHIPPED _____ **DATE REQ'D** _____

*Case will ship to arrive before the APPOINTMENT DATE.
 Please account for transit time to/from your office.

APPT DATE* _____

APPROVED TO CHARGE FOR EXPEDITED SHIPPING TO RETURN ON DATE NEEDED

CALL ME REGARDING CASE EMAIL ME REGARDING CASE

CUSTOM DESIGNED ANCHORAGE

Sintered Bands
 Sintered Crowns - Full Coverage
 Sintered Crowns - Removed Occlusal


MARK ALL TEETH TO BE SEATED

R	7	6	5/e	4/d	d/4	e/5	6	7	L
	7	6	5/e	4/d	d/4	e/5	6	7	

CUSTOM ANCHORAGE DESIGN
 Custom anchorage can either be done in single tooth options (standard on the 6's) or in groups, such as a grouping of the bicuspid to the molars. There is no additional cost by customizing the grouping of your anchorage.

ATTACHMENT OPTIONS

Seating Lugs
 Weld Brackets to Bicuspid Bands
 .018 Buccal Tubes .022 Buccal Tubes



SPECIAL INSTRUCTIONS

DR. SIGNATURE: _____

LICENSE NUMBER: _____ **EXPIRES:** _____